




For Accurate & Trustworthy Blood Glucose Results

# GLUCOCARD™ Σ BLOOD GLUCOSE METER



# SAFE EJECT



Blood Glucose testing now becomes much safer with the 'Safe Eject' Button. No need to handle a used strip. After the test just push the button behind the meter and the used strip drops automatically.



### Biosensor Technology

Eliminates blood contact with the meter. No hassels of disinfecting the meter after each test.



### In-built Thermosensor

Get accurate readings whether performing the test in ICU (air-conditioned environment) or in regular wards (non air-conditioned environment). Glucocard  $\Sigma$  is equipped with a thermosensor, which senses the environment temperature and adjusts to provide temperature-corrected readings.



### 7 Sec Test Time

Saves time when no.of patients to be tested are more.



### Auto-coding

No hassles of remembering to manually change the code every time a new vial is opened or a meter is carried to another ward. Also, no chance of getting in-accurate readings due to a "wrong" code.



### GOD (Glucose Oxidase) Based Enzyme

Ensures only true glucose is measured and other sugar like Dextrose, Fructose / Maltose do not interfere with blood glucose readings.



### 0.5µL Blood Sample


Makes it easier to obtain blood sample from an old patient/ pediatric patient. Now get accurate readings even with a smaller blood drop.

## GLUCOCARD $\Sigma$ Test Strips

**Sensors With Ruthenium Mediator:** Ensures accurate readings even at higher temperatures

**Longer Strip:** Facilitates ease of handling

**Short-fill Detection (Sufficient-sample Indicator):** Eliminates errors or inaccurate readings due to insufficient blood sample.



Use GLUCOCARD Sigma Test Strips with TRUSTLET Safety Lancets for a safer testing experience.



28G x 1.8mm  
28 Gauge thickness

### Safety

Needle is fully shielded before and after use  
Sterilized by gamma irradiation

### Comfortable

High quality tri-bevel with high-speed penetration  
Tri-bevel needle tip keeps good blood flow  
The optimum ergonomics design

### Convenient

Pre-loaded and pressure activated device



# TRUSTLET

Safety Lancet

# 100 Nos.

Safe and convenient single use lancets

For more information contact:

+91 777705331 | customersupport@arkray.co.in

ARKRAY Healthcare Pvt. Ltd. 701-702, Opulence, 6th Road, TPS III, Santacruz (E), Mumbai - 400055, Maharashtra, India.