



- 👍 BROAD RANGE OF HIGH QUALITY PRODUCTS WITH MODERN DESIGN
- 👍 UNCOMPROMISING PERFORMANCE

*Genuine quality speaks for itself*

# ELECTRA BL

## Lab Centrifuge

A STURDY WORKHORSE  
FOR YOUR LABORATORY



- ☒ Maximum Speed: 4000 rpm
- ☒ Motor: DC brushless motor
- ☒ Rotor capacity: 12/18 positions
- ☒ Rotor type: Fixed angle rotor
- ☒ Time range: 0-99 minutes
- ☒ Noise: < 70 db

# €LECTRA BL Swing Out-24T

## MICROPROCESSOR CONTROLLED SWING HEAD CENTRIFUGE



24 TUBES



- ☒ Maximum Speed: 4000 rpm
- ☒ Motor: DC brushless motor
- ☒ Rotor capacity: 24 positions
- ☒ Tube capacity: 5 ml tubes
- ☒ Rotor type: Swing head rotor
- ☒ LCD type preset & time upto 60 min

# LECTRA

Digital Colorimeter  
(AUTO ZERO)

The Pioneer....  
Still the leader....




**NOW TRUELY DIGITAL**  
with modern ergonomic  
& all touch pad.

- ☑ Light: Long LED.
- ☑ Display: Digital display.
- ☑ Filters : 9 filters (410,430,470, 490,520,540,580,610,640)nm.
- ☑ Zero setting: Microprocessor. based auto zero setting.
- ☑ Type: with battery/without battery.

# €LECTRA Blood Mixer

WORLD CLASS QUALITY  
FOR YOUR LABORATORY



- ❑ Rolling speed: 10–30 rpm.
- ❑ Rolling frequency: 5 cycles per min. at 35° angle.
- ❑ MS coated exterior body.
- ❑ PVD coated rolling rods, each of 230 mm.

# €LECTRA

---

## Digital Incubator

A STURDY INCUBATOR  
FOR YOUR LABORATORY




- ☑ Tube positions: 12 tubes
- ☑ Type of tube: Round
- ☑ Temperature: 37°C
- ☑ Power Adaptor : 100 - 200v, 50/60 Hz

# TopPette Micro Pipettes

## PREMIUM QUALITY IMPORTED MICROPIPETTE



- ❑ Premium quality micropipette for high-end laboratories.
- ❑ Accurate, reliable & comfortable for extensive pipetting.
- ❑ Light weight & ergonomic design for user comfort.
- ❑ **Volumes**  
Fixed Volume : 5 $\mu$ l, 10 $\mu$ l, 20 $\mu$ l, 50 $\mu$ l, 100 $\mu$ l, 500 $\mu$ l, 1000 $\mu$ l  
Variable Volume : 5 – 50 $\mu$ l, 2 – 20 $\mu$ l, 10 – 100 $\mu$ l, 20 – 200 $\mu$ l, 100 – 1000 $\mu$ l



For more details:

 [www.peerlessbio.com](http://www.peerlessbio.com)

 [info@peerlessbio.com](mailto:info@peerlessbio.com)

 044 2496 0541

 **TOWARDS**  
**EXCELLENCE** 